

# Hemo-Nate®

Microaggregate Blood Filtration System

Reduce the risk of pulmonary distress and overburdening of underdeveloped neonatal vasculature

Removes small volumes of particulate debris which causes significant changes in pulmonary function of fragile neonates


18 micron filter has high microaggregate retention, yet allows viable blood product to pass through with ease

Complete needleless set available

## INDICATED FOR:

- Neonatal Whole Blood Infusion
- Neonatal Packed Cell Infusion
- Platelet Infusion
- FFP (Fresh Frozen Plasma)
- AHF Concentrates
- Cryoprecipitates
- I.V. Fat Emulsions
- TPN Preparation
- Other Blood Fractionates
- Drugs (as a particulate filter)

\*Not indicated as an anti-bacterial or anti-microbial filter



**Stainless Steel Filter Media** > The electrostatic properties provide exceptional flow, yet is more gentle on blood components than plastic strand fibers.



UTAH MEDICAL  
PRODUCTS INC.



# PRODUCT SPECIFICATIONS

PRODUCT  
NUMBER

DESCRIPTION

PACKAGING

4020005

Complete Needleless Hemo-Nate® Filter Set

50/Box



4020009

Hemo-Nate® Filter

50/Box



4020007

Hemo-Nate® Filter Set

50/Box



4020509

Hemo-Tap®

30/Box



4020609

Hemo-Tap® Valve

30/Box

